

Breezy Spot™ Specifications



Breezy Spot™ extends the capabilities of Breezy Blue™ by automatically recording the exact location of each disinfection cycle—providing hospitals with verified data to support safety and compliance. These small wireless beacons integrate seamlessly with

Breezy Blue and Breezy Cloud™, creating a verified record of disinfection activity across hospital rooms, operating suites, and other critical spaces. Infection preventionists and EVS teams gain auditable, location-based data that simplifies compliance reporting, strengthens accreditation readiness, and builds confidence in cleaning workflows. With Breezy Spot, hospitals can track and document terminal cleanings with ease—reducing manual record-keeping while strengthening accountability and patient safety.



Item	Specification
SKU	BSPT01
Wireless protocol	Bluetooth® LE 5.3
Size	Diameter: 2.19 in. (55.6 mm), Width: 0.91 in (23.1 mm)
Weight	1.44 oz. (40.9 g)
Environmental protection	IP67
Battery	ER14250 (X2), replaceable
Battery Life	Approximately 6 years
Mounting method	3M VHB tape
Operating temperature range	-4° F to 140° F
Operating humidity range	0 to 95% relative humidity, non-condensing
Breezy Spot location - operating distance	Within 7 feet of Breezy Smart Controller during fogging
Breezy Spot location - separation	Greater than 25 feet from other Breezy Spots
Required equipment for location tracking	Breezy Smart Controller (Gen 2), Breezy Blue, Breezy Buggy